

NoteBookMaker®

USER MANUAL v 10.0




Table of Contents

| | |
|---|--------------|
| Welcome and about the NoteBookMaker® | 2 |
| Video Tutorial Support | 3 |
| Installation..... | 4 |
| Using the Software | 5-20 |
| User Set-UP | 5-6 |
| User Login..... | 7-8 |
| The Cover..... | 9 |
| The Page (Text only) no graphics | 10 |
| The Page (Text and Graphics) | |
| Page Headers | 12 |
| Value Lists | 12 |
| Scripts and automation | 12 |
| Search Layout..... | 13 |
| Table of Contents (TOC) | 14 |
| Export – Import- Formats..... | 14 |
| The Librarian..... | 15-17 |
| Card view | 15 |
| Scripts and automation | 16 |
| Table view Card(s) & Book(s)..... | 17 |
| The Security Log | 18 |
| Audit Reports..... | 18 |
| Pasting a Table | 19-20 |
| License | 21-22 |
| Single User | 21-22 |
| Multi-User | 21-22 |

Welcome and thank you for selecting
NoteBookMaker® - The World Leader in Electronic Laboratory NoteBooks.

NoteBookMaker™ is a legal notebook system for laboratory professionals. The solution was produced with FileMaker® Developer database software. It is simple to use and manage. Unlike other database solutions, this product delivers true “what you see is what you get”, a key feature for legal defense and support for 21CFR Part 11 compliance.

NoteBookMaker™ can be setup in single user, multi-user (pier to pier, or client/server). The multi-user mode requires the purchase of FileMaker Pro 7.0 or greater and is best hosted by FileMaker Pro Server. Minimal IT knowledge is required to get started.

It contains 2 main database files a NoteBook, and a Librarian file. The NoteBook file is where the scientist enters their data. The Librarian is where the IT department or administrator sets up users via Library Cards. The NoteBook file looks like a book but has the power of a modern database engine. One can easily create, search, and print scientific data. It uses 2 layouts: Layout 1 text only, and Layout 2 for text and graphics (pictures, sounds and movies).

The direction tool bar controls navigation. Use of the bar controls the correct viewing and printing of data. It was necessary to develop this direction tool bar in order to properly select the correct layout. The record browser is a good visual tool to guide the user through the notebook and list of Cards in the Library



Printing will automatically select the appropriate layout. Data that is ready for patent submission can easily be imaged to pdf, or as a stand-alone data file (similar to this solution) that requires no additional software. This is critical for 21CFR Part 11 compliance.

The default font is Times New Roman - True Type (12 point). It was selected for its legibility both on screen and on paper. The font face is used in all fields accept the graphics field on layout 2. The graphics field can take a paste from any clipboard compatible program. Sophisticated scans, tables, sound, and movies can be placed. Both layouts have a text field for scientific data. Because this solution works like paper and not just a database, it is limited to how many lines of data that can be displayed in an 8.5 x 11 format. In order to provide the legal aspects of an electronic notebook, some pull down menus were stripped of their normal FileMaker Pro functionality. Button commands were substituted for the pull down menus to provide specific access & control for the end user and the IT department. Help files are located on the tabs of each layout:



Tool Tips are available for each button and field label. They can be turned off after the user is familiar with the software.

Video Tutorial Support

This manual is supported by QuickTime videos that can be found by going to <http://notebookmaker.com>. We recommend watching the NoteBookMaker® tutorial videos to fill in any questions that you may have. The videos can easily be included into a company's standard operating procedure training requirements. A list of vides can be found on the next page. More videos will be added over time, so visit the site periodically.

NoteBookMaker for PC and Mac
The World Leader in Virtual NoteBooks

NEWS

- [Free Download](#)
- [NoteBookMaker 8.5](#)
- [Universal Binary Compatible](#)
- [WinXP-Vista](#)
- [MacOSX](#)

FileMaker Solutions Alliance ASSOCIATE

Apple Authorized Business Agent

NoteBookMaker® is a totally secure and completely Virtual – Electronic Laboratory Notebook. The software replaces your traditional paper notebook with a modern system that protects your intellectual property and process. It even has a Virtual Librarian (VL) to take care of NoteBook Administration. Simply, copy-paste your data from any program on to a page. Use one notebook per author or many. You can even merge multiple notebooks into one notebook. You have total flexibility!

NoteBookMaker® has all the benefits of a database such as: searching, printing, emailing, importing and exporting, MS Office and Adobe Acrobat integration. And best of all, your data is nicely presented in a 8.5"x11" notebook format. The software is used by scientists, engineers & many others who need a notebook system that has true document control. Simply **Acquire** your Personal Library Card from your System Administrator (your Librarian), **Enter** a User Name & Password in the Security Log, & **Start** generating legally defensible, 21 CFR Part 11 Compliant data. No IT support required! Low system overhead! We did our job so you can do yours!


Personal Library Card → **Security Log** → **NoteBookMaker™**

NoteBookMaker® has several advantages over paper systems. Firstly, the ability to perform a search on one or many parameters to generate an intelligent found set. With this ability the scientist is able to correlate their findings with other colleges or past work. Gone are the days of wasting company resources due to recreating entire.


Installation The program is easy to set up and use. We kept it simple so that it would work as designed, like a laboratory notebook. So let's get started.

The installation program **NoteBookMakerInstaller** will prompt you to install the application and to read the **license agreement**. Installation will default to the program folder directory on the hard drive, and will also place an icon on your desktop. On a Mac, look inside your applications folder and drag it to your dock for easy access.

The installation process will install the NoteBookMaker 10.x application, 10 NoteBooks, a Training-NoteBook, & a Template-NoteBook. On the PC, there are an additional 30 minor files, and on the Mac there is an extension folder. All the files are in the same directory of the program folder. See the section at the end for a complete listing of all files installed. Because only one directory is installed, backup & archiving is a snap. See the sections on backup & archiving strategies.




Directory view of the NoteBookMaker installation. Mac version is shown.



A list of Training videos demonstrating best practice approaches.

Videos require QuickTime Pro Player.

Learn how to create a Library Card, Edit a card, and manage users and notebooks.



A quick tour of the notebook cover and page.

Using the Software

Before we can begin entering data into our 1-Notebook, we **must** setup the Librarian with a user and a witness. This is a formal process for each user.

Note: If you are using the demonstration software then you can skip this section for now. Simply open the 1-NoteBookMaker.fp7 file and when prompted, type in **Albert Einstein** for the User and **emc2** for the password. This will log you in as Albert Einstein. Later you can go to the Card View in the Library and create your own card. When you are ready to purchase a full license then follow the steps below to gain access to the Librarian.

User Set-Up (The Librarian)

Typically, the Librarian is a person who manages the system but they may also be a user as an author, witness or validator. In the single user version, the user is also the Librarian. The Librarian will administer the notebooks and Library Cards.

First, we must open the Librarian data file. If this is your first attempt, you will see 9 Cards for the Training-NoteBook. You can use their ID to explore the functions of the Notebook system.

Simply click the corresponding Notebook ID,


The screenshot shows the 'Librarian' application window. At the top, there's a toolbar with icons for navigation and a 'Records' section showing '9 Total (Sorted)'. Below this is a 'Layout' dropdown set to 'Book List' and a 'View As' section with icons for list, grid, and preview views. The main area contains a table with columns: View, CARD, CARDS, BOOKS, HELP, CLOSE, Edit, Create or Edit Cards, Demo Mode, Search, and FIND. The table lists 9 training notebooks, each with a 'Click' icon, 'Note Book ID', 'Author Name', 'Picture', 'Card #', 'Security', 'Status', 'Description', and 'Department'.

| Click | Note Book ID | Author Name | Picture | Card # | Security | Status | Description | Department |
|-------|-------------------|-------------------|---------|--------|----------|--------|-------------|------------|
| > | Training-NoteBook | Albert Einstein | | 1 | 3 | Active | Research | Physics |
| > | Training-NoteBook | Charles Darwin | | 2 | 3 | Active | Research | Biology |
| > | Training-NoteBook | Niels Bohr | | 3 | 3 | Active | Research | Chemistry |
| > | Training-NoteBook | James Watson | | 4 | 3 | Active | Research | Biochemi |
| > | Training-NoteBook | Francis Crick | | 5 | 3 | Active | Research | Biochemi |
| > | Training-NoteBook | Linus Pauling | | 6 | 3 | Active | Research | Chemistry |
| > | Training-NoteBook | Maurice Wilkins | | 7 | 3 | Active | Research | Quality |
| > | Training-NoteBook | Rosalind Franklin | | 8 | 3 | Active | Research | Molecular |
| > | Training-NoteBook | Marie Curie | | 9 | 3 | Active | Research | Physics |

and when prompted, type in the user name and password, such as Albert Einstein, emc2.


Now that you have familiarized yourself with the NoteBook, you may want to create a Card for your specific ID. To do this you will need to access the Card layout.

If you have made cards before and want to edit or create new ones, then click the “Create or Edit Cards” button to gain access to the Card layout:



The system comes with the Librarian in Demo Mode, and thus you will not need to type in the administrative user name and password. However, once the Librarian has been registered, then you will need to type in the Librarian user name and password. The registration license is sent via email at time of purchase. Registering the Librarian secures access to the Card layout where all passwords are stored.

For the sake of discussion, let us assume you are in Demo Mode and you have access to the Card layout so you can administer the card system and get familiar with managing the users for the author, witness and validator roles. Once you click the Create or Edit Cards button, a Registration dialog is presented. Select Demo as shown below.



Once clicked, the Card layout is displayed as shown on the next page.

Each system comes with 9 cards that can be edited or deleted. In this case, we will click the New Card button or just modify an existing card.

Librarian

View CARD CARDS BOOKS HELP CLOSE Edit NEW CARD DUPE CARD PHOTO Search

Find FIND-ALL SORT OPEN BOOK

Secure Re-Login PASSWORD SHARING SECURITY

host Privileges

NoteBookMaker® V8.5 - 2007

Mandatory

LIBRARY CARD # 1 # Status

Note Book ID 1-NoteBook 1 Active

Author Name Your Name Here

Password Your Password Here

Security Clearance 3 1=Author, 2=Witness, 3=Author/Witness
6=Browse, 10=Author/Witness/Validator

Department Biochemistry

Company Name Your Organization Here.

Picture

Add your Photo ID

Non-Mandatory

Employee ID Optional information

Group Leader

Issued By Name Stephen Arpie

Issued Date, Time 12/17/2006 10:10:23 PM

Returned To

Returned Date, Time

Modified Date, Time 12/17/2006 10:11:31 PM

Comments

Add your comments here:

Completely fill out the card for the first user. Include a photo to yield a visual ID of the user. Use a security setting of 3 to allow the new user to witness the first page of the 1-NoteBook. Create a second card for 1-NoteBook that will be used for witnessing. This will be for the person that will review the notebook so here too you will need a security setting of 3. Create as many cards as you have users. It is a good idea to have all your photos and user information before beginning this process. The “host” **Password** should be changed after the system is set-up. See the section on Passwords and Security. Once you have made a Library Card, the user can open a book and login.

In multi-user mode, the Librarian administrator will keep the Librarian.fp7 file open for users to access during their normal login process. It is highly recommended that a dedicated server host the Librarian, and all notebook files. The server should not be used for routine data entry of notebooks, but merely hosting and system administration.


User Login

Now that you have your Card completed, you can open the notebook that it references.

To open the 1-NoteBook, click the  button or the Open Book button.



Alternatively, you click out to the Books layout and select as show page 5. In either event, you will be prompted by this dialog:



Select User. User access allows for adding pages and content to the notebook system.


The next step will ask for the user name and password that was typed on the card.

When your 1-NoteBook.fp7 file boots up, you will be prompted for:

User Name “your name”

Password “your password”

This is the information that you typed into the Librarian when you set up your card.


A screenshot of a 'Security Log In' dialog box. The dialog has a blue title bar with the text 'Security Log In'. Below the title bar, the text 'User Name and Password' is displayed. There are two input fields: the first is labeled 'User Name' and contains the text 'Your Name Here'; the second is labeled 'Password' and contains a series of asterisks '*****'. At the bottom of the dialog, there are two buttons: 'Cancel' and 'OK'.

During the login process, the Librarian Security table records the information that the user types. The information typed is compared to the information that was set-up in the Librarian by the administrator/host. If the information is a complete match, then the user is presented with their notebook data file (1-NoteBook.fp7). If the information is not a match, then the user is rejected and has a total of 3 attempts before all the data files are closed.

Upon successful login, the User will appear on the cover as you are presented with your 1-NoteBook (See next section). Their first assignment should be to witness the initial page. This tells management that they have read and understood that this is a secure electronic notebook system.

Click around to get the feel of the application. The user can now start creating and posting pages. On a routine schedule, the person designated as the witness for this 1-NoteBook should review and post their signature to completed pages. The witness will need a Library Card with a security level of 2, 3 or 10.

The Cover



The cover displays the photo identification of the current log in user. In upper right hand corner is a control menu that performs the following functions:



Open will open the notebook to the last page.

TOC opens to the table of contents view.

Library opens the Librarian file.



Help brings up the help information for the cover layout.

Close will close the notebook but keep the application running.

Exit will close all active books and quit the application.

The Page: (TEXT only) no graphics


When you create a new page, you will be prompted to choose between 2 layout styles.



From here you will perform the functions that are listed on the button pallets, such as **Copy**, **Paste**, **Print**, **Go To** navigation, & most importantly, issue **Legal** functions such as **New Page**, **Edit Page**, **Post Author**, **Post Witness**. The security clearance that was set up in the librarian will determine what legal functions that are permitted. There are other constraints as well, such as: an author cannot witness their own page; a witness cannot witness a page that has not been authorized; a “host” cannot authorize or witness any pages, but merely browse. These constraints are essential to legal notebook system. Detailed information can be found in the Help Tabs of each palette.

The Page: (Text and Graphics)

The first layout only involves a body that has one big text area. This text only layout is typically used for free form typing and pasting in text from other programs. If one wants to preserve the




formatting or paste in pictures, then Layout 2 “Text & Graphics” is required. The graphics field above can support many popular formats and thus most popular scientific and graphic programs. Supported formats include: copy/paste, .tiff, .eps, .pdf, .jpg, .gif, and many others. Sound (.aiff, .mp3) and Movie formats are also supported.

Page Headers

As with any page, you should fill out the header fields so you can do proper searches. For example, the Project ID field can be loaded with default values so all your users conform to proper language. Attachment can hold file archives of any type of data.

| | | | |
|---|--------------|----------------|----------|
| Company-Dept Chemicals - R Us, Inc. 1313 Mockingbird lane - Physics | | | |
| Project ID | Double Helix | Page # | 17 |
| Project Name | Heredity | Reference | 16 |
| Key Word | DNA | Type | Original |
| Attachment | | Correction Pg# | |

Value List can be Exported and modified in a spreadsheet program and then copied back into any




notebook. The last item Export Value Lists is what is used to export not only Project ID, but also Project Name, and Compound Name.

The files created can be found in the parent folder of your 1-NoteBook data file.:

Project Name Value List.tab
Page Type Value List.tab
Project ID Value List.tab

Search Layout

With the header fields in combination with the footer and main body text fields, a user can perform very detailed searches.



The buttons to the right, under **Search** tab: **Show All**, **Find**, **Re-Find**, and others are convenient ways to conduct searches. Navigational controls are listed under the **Go To** tab, while the **Edit** functions are listed at the top left.

Enter your search criteria, create multiple requests, expand, and constrain to get the page records that are of interest.


After you have completed your search, you can print, or report on your found set.

It is important to note that only one book can be searched at one time. If the user wishes to search multiple books then they have two options. The first option is to use merge existing notebooks into one notebook. This can be accomplished by importing multiple books into an empty Template-NoteBook.fp7. You can rename the file before or after the process as long as you have a library card that grants access. Once all the file(s) are imported, then a search will be for all the data. Alternatively, you can export each individual notebook as an Adobe Acrobat Portable Document “.pdf”. The network version can export directly to this format. Later, you can import each pdf into a master pdf and use Acrobat’s search routines.

Either way, exporting and importing into a master notebook or master pdf are great ways to search on all your organizations data. And, you do not need an IT department to accomplish it. All these functions are straightforward, end user skills.

Table of Contents (TOC)

Each notebook will automatically build a table of contents. The **Go To** tab provides a great way to perform quick reviews of found data sets. Each field can be sorted to



| Project ID | Project Name | Key Word | Page # | Created Date |
|---------------------------------------|------------------------------------|--|--------|-------------------------|
| Author | Posted Date | Witness | File | |
| ADH James Watson | ADH 12/4/2006 | ADH Francis Crick | 8 | 12/4/2006 1-NoteBook |
| Alcohol dehydrogenase James Watson | Alcohol dehydrogenase 12/4/2006 | Alcohol dehydrogenase Francis Crick | 10 | 12/4/2006 1-NoteBook |
| Catalase James Watson | Catalase 12/4/2006 | Catalase Francis Crick | 6 | 12/4/2006 1-NoteBook |

generate a specific order. The fields displayed are: Project ID, Project Name, Compound Name, Page #, File Name, Author, Posted Date, Witness and Date Created. The reports can be output to any device or electronic format.

An advanced technique involves combining several notebooks into one. The text fields can be exported from each notebook file and imported into a master notebook. This way, all fields can be searched across all notebooks. The TOC will display a master page directory for all your entries.

The export can be .TAB, .XML, or .fp7 formats. The .fp7 is the preferred way to export and import for this feature. It will retain the original formatting and file names.



The Librarian (Detail)

As mentioned in the Setup pages, the Librarian is the place to register your users. Each User shall receive a card with the mandatory fields listed below. The **Go To** tab shall display the following information:

Librarian

Layout: Card View

Record: 1

Total: 1

Unsorted

Mandatory

LIBRARY CARD # 1 # Status

☒ Note Book ID 1-NoteBook 1 Active

Author Name Your Name Here

Password Your Password Here

Security Clearance 3 1=Author, 2=Witness, 3=Author/Witness
6=Browse, 10=Author/Witness/Validator

Department Biochemistry

Company Name Your Company Name Here

Picture

Add your Photo ID

Non-Mandatory

Employee ID Optional Information

Group Leader

Issued By Name Stephen Arpie

Issued Date, Time 12/18/2006 9:49:06 AM

Returned To

Returned Date, Time

Modified Date, Time 12/18/2006 9:49:47 AM

Comments

Add your comments here:

Secure

Re-Login


PASSWORD

SHARING

SECURITY

host Privileges

NoteBookMaker® V8.5 - 2007



Card View button brings the Librarian Administrator to the card layout. Create and edit cards here. This layout requires Librarian Administration privileges. Normal User privileges will not allow access to the layout. This is because User Names and Passwords are kept here.

Cards View button brings the Librarian Administrator to the table view of found set of cards. Similar to the card view layout, no access is allowed for normal users. This layout is a great way to view all or some of the cards in the Library.

Book View button brings the user to the book view. This is the entry point to the Library. From here you can choose a book to open or Re-Login as the Librarian as the user.

Help button opens the general help layout.

Close will only close the Librarian data file, and not other active files. This will prevent any user from login into any notebook. The Librarian must be open for others to access their notebook.



New Card button creates a new blank card that needs to be filled out before using.

Dupe Card button duplicates the existing library card.

Photo allows the user to directory locate a file that is to be imported in to the user picture. The file format is compatible bmp, cgm, drw, gif, macpaint, pcx, pic, jpg, lotus SLD, tiff, metafile, eps, and epsf. Be sure to select the document preference for Store Compatible Graphics before importing each graphic into your notebook. The preference file will direct the notebook to store cross-platform graphics.



Find places the user in the search mode. Additional criteria commands are available to assist the user in completed a search.

Find-All will show all the card entries in the database.

Sort button allows the user to sort by specified criteria

Open Book button opens the notebook listed in the Note Book ID field. It is designed to provide a fast open routine for systems that contain a large number of notebooks.



button will also open the listed NoteBook.




Re-Login allows the Librarian to log out “leave” the system in a secure mode. Securely type at the prompts Account Name and Password.

Password button provides the identification signature in conjunction with the user’s name.

Sharing allows the Host to change the availability of the Librarian to 1 of 3 states: 1. Single User, 2, Multi-User, and 3, Multi-User hidden. Multi-User is the state necessary for networking.

Security button will open the security log layout and allow the Librarian administrator to view the audit trail report.

In **Card View** mode, the found set of records is **navigated** by clicking on the record index in the upper left hand corner of the screen.



Some additional features are the ability to sort your card set by clicking on the field headers such as: Library Card #, Note Book ID. A complete listing of automation scripts is shown here to the right above.



button will also open the listed NoteBook.

Table View Mode in Detail (Card List for Librarian)

It is obvious that The Librarian is the place to hold all your user information, but there are some other interesting features that system administrators shall appreciate.

| Card # | Author Name | Picture | Note Book ID | # | Pass Word | Security | Status | Department | Company Name | Issued Date |
|--------|-------------------|---------|--------------|---|-----------|----------|--------|-------------------|------------------------|-------------|
| 1 | Albert Einstein | | 1-NoteBook | 1 | emc2 | 10 | Active | Physics | Chemicals - R Us, Inc. | 3/3/200 |
| 2 | Charles Darwin | | 1-NoteBook | 1 | emc2 | 10 | Active | Biology | Chemicals - R Us, Inc. | 3/3/200 |
| 3 | Niels Bohr | | 1-NoteBook | 1 | emc3 | 10 | Active | Chemistry | Chemicals - R Us, Inc. | 4/22/200 |
| 4 | James Watson | | 1-NoteBook | 1 | emc2 | 10 | Active | Biochemistry | Chemicals - R Us, Inc. | 4/22/200 |
| 5 | Francis Crick | | 1-NoteBook | 1 | emc2 | 10 | Active | Biochemistry | Chemicals - R Us, Inc. | 4/22/200 |
| 6 | Linus Pauling | | 1-NoteBook | 1 | emc2 | 10 | Active | Chemistry | Chemicals - R Us, Inc. | 5/9/200 |
| 7 | Maurice Wilkins | | 1-NoteBook | 1 | emc2 | 10 | Active | Quality | Chemicals - R Us, Inc. | 12/1/200 |
| 8 | Rosalind Franklin | | 1-NoteBook | 1 | emc2 | 10 | Active | Molecular Biology | Chemicals - R Us, Inc. | 1/20/200 |
| 9 | Marie Curie | | 1-NoteBook | 1 | emc2 | 10 | Active | Physics | Chemicals - R Us, Inc. | 2/17/200 |
| 10 | Stephen J. Arpie | | 1-NoteBook | 1 | emc2 | 10 | Active | Quality | Chemicals - R Us, Inc. | 11/12/200 |

A very useful tool built into The Librarian is the **Cards View** mode. It is different than the **Card View** mode primarily because it shows you all the users in a found set. If we issue the **Find All** command before going to table view we will see a complete listing. The Table view can be sorted by clicking on the individual headers, such as Author Name or NoteBook ID. This is extremely helpful when you need to identify who and what someone has access to use or has used. This is a major component of 21CFR Part 11 compliance. System administrators can quickly modify access privileges on the found set by applying replace content command on a particular field. Care must be given before issuing the function. Modifications can't be undone unless they are manually type set.

Book List View Mode in Detail (Book List for users)

Books List mode is the main entry point in the Library.

| Note Book ID | # | Author Name | Picture | Card # | Security | Status | Department | Issued Date |
|--------------|---|-------------------|---------|--------|----------|--------|-------------------|-------------|
| 1-NoteBook | 1 | Albert Einstein | | 1 | 10 | Active | Physics | 3/3/200 |
| 1-NoteBook | 1 | Charles Darwin | | 2 | 10 | Active | Biology | 3/3/200 |
| 1-NoteBook | 1 | Niels Bohr | | 3 | 10 | Active | Chemistry | 4/22/200 |
| 1-NoteBook | 1 | James Watson | | 4 | 10 | Active | Biochemistry | 4/22/200 |
| 1-NoteBook | 1 | Francis Crick | | 5 | 10 | Active | Biochemistry | 4/22/200 |
| 1-NoteBook | 1 | Linus Pauling | | 6 | 10 | Active | Chemistry | 5/9/200 |
| 1-NoteBook | 1 | Maurice Wilkins | | 7 | 10 | Active | Quality | 12/1/200 |
| 1-NoteBook | 1 | Rosalind Franklin | | 8 | 10 | Active | Molecular Biology | 1/20/200 |
| 1-NoteBook | 1 | Marie Curie | | 9 | 10 | Active | Physics | 2/17/200 |
| 1-NoteBook | 1 | Stephen J. Arpie | | 10 | 10 | Active | Quality | 11/12/200 |

Here the users can easily find there notebooks and open them. No passwords are viewable from this layout. The Librarian in conjunction with the Security Log


provide full security & audit reporting which are key to 21CFR Part 11.

With these tools, you will easily satisfy your auditor's requirements.

The Security Log

From the Card view of the Librarian, the Administrator can readily obtain information about who has opened a notebook and when access has occurred. Also, if there are any failed attempts at login, then they too are recorded by the Security Log. To view this information one simply clicks **Security** button from the Library Card View – Secure scripts palette. Only the Librarian administrator should view this information because it contains all the passwords entered by each user.

The system administrator is shown the cabinet view:




There is a control palette with an **Audit** button that the administrator/host clicks on to gain access to the report.

By clicking the **Find-All** button, a complete report of all activity is generated. Modifying

Dates, names, passwords, and sorting too can customize the report. The audit trail is another key requirement of 21 CFR Part 11 compliance.

During the login process, the Security database opens along with the Librarian and records the information that the user types. The information typed is compared to the information that was set-

up in the Librarian by the administrator/host. If the information is a complete match, the



user is presented with their notebook data file (1-NoteBook.fp7 in stand-alone mode, and 1-NoteBook.fm7 in multi-user mode). If the information is not a match, then the user is rejected and has a total of 3 attempts before all the data files are closed.


Pasting a Table as Text and as a Graphic.

NoteBookMaker supports Microsoft Excel, Word and PowerPoint files via the clipboard for easy copy-paste functionality.

The advantage of pasting in the table as a graphic element is that the layout is maintained and cannot be changed. This too can be copied out and re-used if required. Any size table can be used in this procedure. Even if the table is very large, it will image out on paper as magazine quality. Alternately, one could copy/paste sections until the complete table is produced. Use the magnification for on-screen viewing of detailed spreadsheets in NoteBookMaker®.


First, you will click the graphic element and copy;

Next, click back on the notebook page, and issue the paste command.



If you want multiple images, such as tables and graphics for the same page, then you will need to shift-click each element, click copy, then bring them over together on the same clip board. They will automatically reduce to fit and center for proper viewing.

This is how the table looks once you paste it into a graphic/text page.



L

icense

[NoteBookMaker® SOFTWARE LICENSE]

PLEASE READ THIS LICENSE CAREFULLY BEFORE USING THIS SOFTWARE. BY USING THIS SOFTWARE, YOU AGREE TO BECOME BOUND BY THE TERMS OF THIS LICENSE. IF YOU DO NOT AGREE TO THE TERMS OF THIS LICENSE, DO NOT USE THIS SOFTWARE AND PROMPTLY RETURN IT TO THE PLACE WHERE YOU OBTAINED IT FOR A FULL REFUND.

The enclosed "NoteBookMaker®" electronic laboratory notebook program ("Software") is licensed, not sold, to you by the author for use only under the terms of this License, and the author reserves any rights not expressly granted to you. You own the media on which the Software is recorded or fixed, but the author and its licensors retain ownership of the Software itself.

1. License. This License allows you to:

(a) Use one copy of the Software on a single computer at a time. To use the Software means that the Software is either loaded in the temporary memory (i.e., RAM) of a computer or installed on the permanent memory of a computer (i.e., hard disk, etc.).

(b) Make one copy of the Software in machine-readable form solely for backup purposes. As an express condition of this License, you must reproduce on each copy any copyright notice or other proprietary notice that is on the original copy supplied by the author.

(c) Stand alone single user license (with demonstration data) can be permanently transferred under this License to another party by providing to such party all copies of the Software licensed under this License together with a copy of this License and the accompanying written materials, provided that the other party reads and agrees to accept the terms and conditions of this License. Transfer of rights only applies to the stand-alone single user version. Transfer of the multi-user license does not apply, (see next section).

(d) The Stand-alone single user license (without demonstration) cannot be permanently transferred under this License to another party. The license requires each unique user to purchase a separate license. Thus, each unique user must be registered via the Librarian to gain access. The license is granted per user/company. If a user changes company employment, they must purchase a new license to use the system in the new company. A license from one company to another is not permitted. Agreements between NoteBookMaker, LLC and the purchaser must be made in advance of sale. NoteBookMaker, LLC retains the right to audit the user company once per year for compliance of the end user license agreement. Contact notebookmaker@mac.com to initiate the number of licenses required to outfit your organization.

(e) The multi-user license requires each unique user to purchase a separate license unless otherwise noted in the sales agreement. Thus, each unique user must be registered via the Librarian to gain access. The license is granted per user/company. If a user changes company employment, they must purchase a new license to use the system in the new company. A license from one company to another is not permitted. Agreements between NoteBookMaker, LLC and the purchaser must be made in advance of sale. NoteBookMaker, LLC retains the right to audit the user company once per year for compliance of the end user license agreement. Contact notebookmaker@mac.com to initiate the number of licenses required to outfit your organization.

2. Restrictions. The Software contains trade secrets in its human perceivable form and, to protect them, you may not REVERSE ENGINEER, DECOMPILE, DISASSEMBLE OR OTHERWISE REDUCE

THE SOFTWARE TO ANY HUMAN PERCEIVABLE FORM. YOU MAY NOT MODIFY, ADAPT, TRANSLATE, RENT, LEASE, LOAN OR CREATE DERIVATIVE WORKS BASED UPON THE SOFTWARE OR ANY PART THEREOF.

3. Termination. This License is effective until terminated. This License will terminate immediately without notice from the author or judicial resolution if you fail to comply with any provision of this License. Upon such termination you must destroy the Software, all accompanying written materials and all copies thereof, and Sections 6, 7 and 8 will survive any termination.

4. Passwords. The Software licensed to you contains features, which are protected by certain passwords. You are not authorized to enter, remove or change such passwords, and only the author is authorized to enter, remove or change such passwords.

5. Export Law Assurances. You agree that neither the Software nor any direct product thereof is being or will be shipped, transferred or re-exported, directly or indirectly, into any country prohibited by the United States Export Administration Act and the regulations there under or will be used for any purpose prohibited by the Act.

6. Limited Warranty. The author warrants for a period of ninety (90) days from your date of purchase that (i) the media on which the Software is recorded will be free from defects in materials and workmanship under normal use, and (ii) the Software as provided by the author will substantially conform to the author's published specifications for the Software. The author's entire liability and your sole and exclusive remedy for any breach of the foregoing limited warranty will be, at the author's option, replacement of the media, refund of the purchase price or repair or replacement of the Software.

THIS LIMITED WARRANTY IS THE ONLY WARRANTY PROVIDED BY THE AUTHOR AND ITS LICENSORS EXPRESSLY DISCLAIM ALL OTHER WARRANTIES, EITHER EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE WITH REGARD TO THE SOFTWARE AND ACCOMPANYING WRITTEN MATERIALS. BECAUSE SOME JURISDICTIONS DO NOT ALLOW THE EXCLUSION OR LIMITATION OF IMPLIED WARRANTIES, THE ABOVE LIMITATION MAY NOT APPLY TO YOU.

7. Limitation of Remedies and Damages. In no event will the author, its parent or subsidiaries or any of the licensors, directors, officers, employees or affiliates of any of the foregoing be liable to you for any consequential, incidental, indirect or special damages whatsoever (including, without limitation, damages for loss of business profits, business interruption, loss of business information and the like), whether foreseeable or unforeseeable, arising out of the use of or inability to use the Software or accompanying written materials, regardless of the basis of the claim and even if the author or the author's representative has been advised of the possibility of such damage. The author's liability to you for direct damages for any cause whatsoever, and regardless of the form of the action, will be limited to the greater of \$500 or the money paid for the Software that caused the damages.

THIS LIMITATION WILL NOT APPLY IN CASE OF PERSONAL INJURY ONLY WHERE AND TO THE EXTENT THAT APPLICABLE LAW REQUIRES SUCH LIABILITY. BECAUSE SOME JURISDICTIONS DO NOT ALLOW THE EXCLUSION OR LIMITATION OF LIABILITY FOR CONSEQUENTIAL OR INCIDENTAL DAMAGES, THE ABOVE LIMITATION MAY NOT APPLY TO YOU.

8. General. This License will be construed under the laws of the State of Connecticut, except for that body of law dealing with conflicts of law. If any provision of this License shall be held by a court of competent jurisdiction to be contrary to law, that provision will be enforced to the maximum extent permissible, and the remaining provisions of this License will remain in full force and effect. If the Software is supplied to the United States Government, the Software is classified as restricted computer software as defined in clause 52.227-19 of the FAR. The United States Governments rights to the Software are as provided in clause 52.227-19 of the FAR.